carefix®

EN: Instructions for Use | **AR:** تعليمات الاستعمال | **BG:** Инструкции за употреба | **CS:** Návod k použití
DA: Brugsanvisning | **DE:** Gebrauchsanweisung | **EL:** Οδηγίες χρήσης | **ES:** Instrucciones de uso | **FI:** Käyttöohjeilla
FR: Notice d'utilisation | **HR:** Upute za uporabu | **HU:** Használati útmutató | **IT:** Istruzioni per l'uso | **JA:** 使用説明
KO: 사용 지침 | **NL:** Gebruiksaanwijzing | **NO:** Bruksanvisning | **PL:** Instrukcja użytkownika | **PT:** Instruções de
utilização | **SK:** Návod na použitie | **SL:** Navodila za uporabo | **SR:** Uputstvo za upotrebu | **SV:** Bruksanvisning
ZH: 使用说明

EN: Post-Surgical Bras | **AR:** حمالات الصدر المخصصة للارتداء بعد العمليات الجراحية
BG: Следоперативен сутиен | **CS:** Pooperační podprsenky | **DA:** Post-operative bh'er | **DE:** Postopera-
tive BHs | **EL:** Μετεγχειρητικοί στηθόδεσμοι | **ES:** Sujetadores postquirúrgicos | **FI:** Rintaliivit
leikkauksen jälkeiseen käyttöön | **FR:** Soutiens-gorge post-chirurgicaux | **HR:** Postoperativni grudnja-
ci | **HU:** Műtét utáni melltartók | **IT:** Reggiseni post-chirurgici | **JA:** 術後用ブラジャー | **KO:** 수술후 보정용 브
라 | **NL:** Postoperatieve beha's | **NO:** Postoperative BHer | **PL:** Biustonosze pooperacyjne | **PT:** Sutiões
pós-operatórios | **SK:** Pooperačné podprsenky | **SL:** Modrčiki za uporabo po kirurškem posegu | **SR:** Post-hirurški
grudnjaci | **SV:** Postoperativa brösthållare | **ZH:** 术后文胸

EN: The device in primary packaging is prepared for direct clinical use | **AR:** مَادِخْتَسَالِا لِحَاجَاتِ مَسَائِلِ اسْتِوَاعِ يَفِ زَاهِلِ الْاَنْزِيَهِيَه جَتِي حَتِي
رَشَابِلِا يَكْتِنِي لِكُلِّا
BG: Устройството в първична опаковка е готово за директна клинична употреба | **CS:** Zdravot-
nický prostředek v primárním balení je připraven k přímému klinickému použití | **DA:** Produkt i primær emballage forberedt til direkte medicinsk brug | **DE:** Produkte in
einer Primäverpackung sind gebrauchsfertig und können direkt für eine klinische Nutzung verwendet werden | **EL:** Η συσκευή στην αρχική συσκευασία
παρασκευάζεται για άμεση κλινική χρήση | **ES:** El producto en su envase primario está destinado al uso clínico directo | **FI:** Kuluttajakauksessa
oleva tuote on valmis suoraan kliniseen käyttöön | **FR:** Dans son emballage primaire, le dispositif est préparé pour une utilisation clinique directe | **HR:** Proizvod u
primarnom pakiranju pripremljen je za izravnu kliničku uporabu | **HU:** Az elsődleges csomagolású eszköz közvetlen klinikai felhasználásra készül | **IT:** Il presidio presente
nella confezione è adatto all'uso clinico diretto | **JA:** 本製品は直接包装されており、臨床にもご使用いただけます。 | **KO:** 기본 포
장재의 기기는 직접 임상적 사용을 위해 준비되어 있습니다 | **NL:** Het hulpmiddel in primaire verpakking wordt voorbereid voor direct klinisch
gebruik | **NO:** Enheten i primærpakningen er klargjort for direkte klinisk bruk | **PL:** Wyrób w opakowaniu podstawowym jest przygotowany do bezpośredniego użytku
klinicznego | **PT:** O dispositivo na embalagem primária está preparado para utilização clínica direta | **SK:** Pomůcka v spotřebiteľskom obale je pripravená na priame
klinické použitie | **SL:** Naprava v osnovni embalaži je pripravljena na neposredno klinično uporabo | **SR:** Uređaj u primarnom pakovanju pripremljen je za direktnu
kliničku upotrebu | **SV:** Produkten i primärförpackningen är förberedd för direkt klinisk användning | **ZH:** 初级包装中的本品可直接用于临床。

3211 Post-Op Bra, Bianca | 3212 Post-Op Bra, Bella | 3270 Post-Op Bra, Anna | 3291 Post-Op Bra, Alice | 3342 Post-Op Bra, Sophia | 3343 Post-Op
Bra, Mary | 3373 Post-Op Bra, Marianne | 3444 Post-Op Bra, Ava | 3471 Post-Op Bra, Lisa | 3609 Post-Op Bra, Bonnie | 3831 Post-Op Bra, Bree
3954 Post-Op Bra, Alma | 3086 Comfort Bra, Catia | 5006 Comfort Bra, Celia 5007 Comfort Bra, Celia

Intended purpose (Post-Op Bras): Fixate and support the breasts and protect wound and scar following breast surgery. Fixate bandages or other medical devices following breast surgery.

Intended purpose (Comfort Bras): Fixate and support the breasts and protect wound and scar when full movement of arm is regained following breast surgery and in radiotherapy period.

Users: Women, following breast surgery as described under Indications.

Indications

For use 24 hours per day following the surgeries:

3211: Mastectomy, lumpectomy, in the radiotherapy period

3471, 3842: Mastectomy, lumpectomy, breast reduction, breast reconstruction

3373, 3343: Mastectomy, lumpectomy, breast reduction, breast reconstruction, breastlift (mastopexy), SN biopsy and axillary dissection of lymph nodes, breast augmentation

3609, 3212, 3291: Mastectomy, breast reduction, breast reconstruction, breastlift (mastopexy), implant removal and replacement surgery (capsulectomy) breast augmentation

3270, 3841, 3954: Breast augmentation, breastlift (mastopexy), mastectomy, lumpectomy, breast reduction, breast reconstruction, implant removal and replacement surgery (capsulectomy)

3444: Breast augmentation, breastlift (mastopexy), breast reduction, mastectomy, breast reconstruction

3086, 5006, 5007: Mastectomy, lumpectomy, in the radiotherapy period

Use the device according to indications to prevent injury or damage.

Contraindications

Relative contraindications:

- Very sensitive scar and skin tissue in the operated area
- Skindiseases
- Radiodermatitis
- Contact dermatitis

Clinical benefits

- Fixation of bandages and wound dressings
- Protection of wound/scar
- Support and compression of the operated breast tissue

Material composition

Raw materials and components used in the device are tested for harmful substances according to STANDARD 100 by OEKO-TEX®, product class II, and relevant Biological Evaluation (ISO 10993). The device and packaging comply with REACH.



Products may contain traces of softening and antistatic agents.

Instructions for use

Depending on medical indication, Post-Op Bras can be combined with 8113 Implant Stabilizer band. 8108 Breast Puffs can be used with pocketed Post-Op Bras.

Fitting the bra and accessories

Post-Op Bras

STEPS	PRODUCT(S)	INSTRUCTION
0.	3954	Easily switch between straight or crossed straps before application by detaching the straps at the front, threading through the opposite top eyelet, and securing with the hook & loop closure. 
1.	All Post-Op Bras	If medical condition allows, insert one arm and pull the bra over the shoulder. Then insert the other arm. If necessary, ask clinical staff for assistance.
2a.	3212 3291 3343 3373 3609 3831	 Pull the under bust band together. Fasten the bottom hook first and work up to the top. Use the widest setting that fits immediately following surgery and the narrowest setting when swelling has subsided.
2b.	3211	 Pull the under bust band together. Work from the bottom to the top, pressing firmly on hook & loop to close.
2c.	3270 3342 3444 3471 3954	 Pull the bra together. Fasten the hooks first, followed by the zipper.
3.	All Post-Op Bras	Position each breast correctly in the cup. Avoid pulling the bra down too low.
4.	All Post-Op Bras	If the bra has adjustable shoulder straps, adjust these so the bra; does not irritate the skin, provides support under the breasts and sides and feels comfortable to wear. 

Warnings and precautions

- Follow the Instructions for Use
- Follow your surgeon's recommendations regarding use and fitting
- If any unexpected discomfort or wound irritation occurs, remove the Post-Surgical Bra and contact your surgeon
- The Post-Surgical Bra is not intended for direct contact with injured skin
- Make sure that bandages, textile breast prosthesis and/or drainage tubes are placed correctly to avoid pressure points and/or obstructed flow
- Before X-ray examination remove Post-Surgical Bra except, 3211, 3086 and 5006/5007
- Wearing an undersized Post-Surgical Bra may cause local pressure pain and disruption of blood flow and/or lymphatic flow
- Wearing an oversized Post-Surgical Bra may reduce the clinical effect

Performance characteristics

General performance characteristics

- Compression and support zones: Integrated and targeted compression and support zones.
- 4-way stretch: The bra adjusts to the operated breast along with the patient's changing breast contour after surgery.
- Seamless: No seams in critical areas (cups, sides, back and under bust).
- Flexible sizing: Bra fits several sizes due to 4-way stretch and size concept.
- Breathable: Knitted material is breathable.

Specific performance characteristics

Reporting obligation

If you have any complaints, please contact your supplier directly. In the events of a serious incident that could lead to a significant deterioration of health or to death, contact the manufacturer and the competent authority of the member state. Serious incidents are defined in article 2 No. 65 of Regulation (EU) 2017/745 (MDR).


Lifetime of the device in use

The lifetime of the medical device is determined by natural wear and tear, even when treated correctly and as recommended. In case of signs of deterioration affecting performance and comfort, the device should be disposed of.

Care instructions

Find care instructions in the bra. Close the bra before washing. Wash with similar colours. Wash in a laundry bag. Use of chlorine based products may decrease durability and function.

Comfort Bras

STEPS	PRODUCT(S)	INSTRUCTION
1.	3086 5006 5007	Pull bra over your head or pull it up like underwear.  
2.	3086 5006 5007	 Adjust the bra and ensure nothing is folded.

Textile Breast Prosthesis

STEPS	PRODUCT(S)	INSTRUCTION
1.	3342 3343 3373 3444 3471	 Remove fiberfill from the device until size fits.
2.		 Place the device in the pocket before bra is fitted.

Mounting drainage bulb

STEPS	PRODUCT(S)	INSTRUCTION
1.	3373	 Place drainage bulb (max 500 ml) in strap and place strap on bra. Press firmly on hook & loop.

Damaged packaging

If the packaging is damaged, do not use the product. E.g. if the packaging is visibly damaged and/or the product itself is visibly compromised or contaminated, the product must be discarded.

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Carefix is a registered trademark of Tytex A/S © 2025.
All rights reserved Tytex A/S, Marsvej 6, DK-7430 Ikast
Tytex A/S has been certified according to ISO 9001, ISO 14001 and ISO 13485
TYTEX A/S • Marsvej 6 • DK-7430 Ikast • +45 9660 4200 • www.tytex.com